




OVARIAN CYST




Ingredients

- Two white rubber balloons
- Liquid honey in a cone-nosed bottle
- Air tank
- Scotch tape or suction cup and crocodile clamp
- Lap trainer




1. Put one balloon inside the honey-filled balloon




2. Inflate the balloon which is inside with air to make a sphere with a diameter of approximately 4 cm



3. Tie off the inner balloon tightly






4. Tie off the outside balloon separately



5. Fasten the crocodile clamp to the suction cup. Attach crocodile clamp to the balloon



6. Secure the model to the bottom of the lap trainer with the bottom towards the student
(Scotch tape can be used in the place of suction cup for anchor)

Cover the trainer with skin